



Cryosurvival and fertility of mouse cumulus oocyte complexes by using a Cryoloop method



Veterinary Research
Scholars Program
University of Missouri

Lyndsy Gieche, Cansu Agca, Liga Wuri, and Yuksel Agca

College of
Veterinary Medicine
University of Missouri

Department of Veterinary Pathobiology, College of Veterinary Medicine, University of Missouri, Columbia, Missouri

Introduction:

- ❖ Embryo and gamete cryopreservation is important to preserving rodent strains used as models for animal and human disease research.
- ❖ Cryopreservation of mammalian oocytes is less successful than that of sperm and embryos.
- ❖ Increased oocyte cryopreservation success would provide more complete genome banking of genetically modified rodent strains.

Objectives:

- ❖ To improve oocyte cryopreservation success by comparing mouse oocyte cumulus complexes (COCs) cryo-survival and post-thaw fertility rates between a vitrification procedure with a Cryoloop versus a 0.25 mL French straw.

Materials and Methods:

- ❖ **Animals and COCs collection:** Outbred CD-1 mice (n=22) were superovulated by an intraperitoneal injection of 10 IU PMSG, and a 10 IU HcG injection 48 hrs later.
- ❖ **Vitrification and thawing procedure:** COCs were exposed to PBS containing 15% ethylene glycol (EG) + 20% fetal calf serum (FCS) for 5 min. The COCs were then transferred into PBS containing 30% EG + 0.5 mol/L sucrose + 20% FCS for 1 min, loaded into Cryoloop or 0.25 mL French straw (8 per method) and immediately plunged into liquid nitrogen (LN₂) for vitrification. The vitrified COCs were thawed by transferring them into PBS containing 0.5 M sucrose and 20% FCS for 5 min at 37°C.
- ❖ **Post-thaw viability evaluation:** COCs were transferred into mouse embryo culture medium (FHM) containing 1 mg/mL hyaluronidase to remove cumulus cells surrounding the oocytes. Numbers of oocytes with intact and damaged oolemma were recorded based on morphological evaluation.
- ❖ **In-vitro fertilization (IVF) and embryo culture (IVC):** An IVF procedure was performed to compare the fertility of COCs vitrified by a Cryoloop method versus freshly collected COCs using fresh epididymal sperm. IVC were performed in KSOM media to determine developmental competence of Cryoloop-vitrified COCs up to blastocyst stage.
- ❖ **Statistical Analysis:** The data was analyzed by the general linear model procedure of SAS. The experiments were repeated at least six times for each experiment.




Figure 3: Mouse cumulus oocyte complex




Figure 4: Nylon loop of Cryoloop over mm portion of a ruler; loop diameter is approximately 1 mm wide



Figure 5: Demonstration of ice crystal formation of a water droplet in a Cryoloop plunged into LN₂




Figure 6: Vitrification achieved with a droplet of cryoprotectant in a Cryoloop plunged into LN₂




Figure 7: Mouse oocytes after vitrification and thawing of COCs with the Cryoloop method; Arrow A shows a dead oocyte, determined by the lack of oolemma integrity and cytoplasmic content; Arrow B shows a viable oocyte with intact oolemma




Figure 8: Mouse oocytes after vitrification and thawing of COCs with the 0.25 mL French straw method; Arrow C shows a viable oocyte with intact oolemma; Arrow D shows a dead oocyte, determined by the lack of oolemma integrity and cytoplasmic content




Figure 9: Mouse 2-cell stage embryos resulting from freshly fertilized COCs using fresh sperm




Figure 10: Mouse 2-cell stage embryos resulting after Cryoloop-vitrification and thawing




Figure 11: Mouse blastocyst stage embryos resulting from freshly fertilized COCs



Figure 12: Mouse blastocyst stage embryos resulting after Cryoloop-vitrification and thawing

Table 1. Post-thaw morphologic integrity of mouse oocytes following vitrification with Cryoloop or 0.25 mL French straws (n=8).

| | No. oocytes vitrified | No. oocytes survived (%) |
|----------------------|-----------------------|--------------------------|
| 0.25 mL French straw | 129 | 25 (19±2.0) ^a |
| Cryoloop | 131 | 99 (75±2.0) ^b |

^{a,b} Values with different letters within a column are significantly different (P<0.05).

Table 2. In-vitro fertilization and pre-implantation embryonic development competence of fresh and Cryoloop-vitrified mouse COCs to blastocyst stage (n=14).

| | No. oocytes | No. oocytes fertilized (%) | No. 2-cell (%) | No. blastocysts (%) |
|-----------|-------------|----------------------------|--------------------------|--------------------------|
| Fresh | 82 | 64 (78±2.5) ^a | 63 (99±2.1) ^a | 58 (90±2.8) ^a |
| Vitrified | 115 | 82 (71±2.2) ^b | 64 (78±1.8) ^b | 46 (56±2.4) ^b |

^{a,b} Values with different letters within a column are significantly different (P<0.05).

Results:

- ❖ Cryosurvival rates of COCs vitrified using 0.25 mL French straws were significantly lower (19%) than those vitrified with the Cryoloop method (75%) (P<0.05).
- ❖ In-vitro fertilization rates of COCs vitrified by the Cryoloop method was lower (71%) than those freshly fertilized COCs (78%) (P<0.05).
- ❖ Embryo development rates following IVF of Cryoloop-vitrified mouse COCs (56%) were lower than those freshly collected COCs (90%) (P<0.05).

Conclusions:

- ❖ Despite some reduction, acceptable in-vitro fertility and embryo development rates can be achieved following cryopreservation of mouse COCs with a Cryoloop vitrification method.
- ❖ This method may be a successful way to preserve mouse strains for genome banking.

Future Studies:

- ❖ Future testing will focus on fetal development competence of the embryos that were derived from Cryoloop-vitrified mouse COCs.

Acknowledgements/Funding:

- ❖ The University of Missouri Mutant Mouse Resource and Research Center (U42OD010918)




Figure 1: Outbred CD-1 mice




Figure 2: Arrow A shows cryoprotectant solution inside a 0.25 mL French straw; Arrow B shows a Cryoloop, approximately 2 cm long, inside the Cryoloop chamber, approximately 4 cm long by 1.5 cm in diameter